Pioneer DreamCare Dr Glucose Non Invasive 先锋一言糖 无创测糖仪

Guide

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www.PioneerIoT.com

- The product name and trademark is reserved to the original manufacturer. Please do not use the product in public without authorization.
- 產品名稱及商標均為原廠所有，未經授權，請勿公開使用

- Dr Glucose is a health information detection device. The detected data cannot be used to replace the diagnosis, counsel, or advices provided by doctors and physicians.
- 本產品為優糖精量健康助理非醫療器材，這些數據不可取代醫師的看診、諮詢、診斷或建議的治療。

English instructions / 中文说明

Foreword/前言	1-01
Product Basics/产品说明	1-02
Appearance/外观说明	1-03
Operation Conditions/工作环境	1-04
Notices/用户注意事项	1-05
IOS (Apple or IPAD) Bluetooth connection settings/IOS(苹果或 IPAD)蓝芽联机设定	1-06
Interface / 操作接口	1-08
History/ 历史纪录	1-19
Self-Calibration/ 自校准模式	1-22
Safety Instructions/ 安全说明	1-24

Thank you for Choosing Dr Glucose!

This manual contains comprehensive information about our product, precautions and many more.

Please read it carefully before using.
We hope you enjoy our product.

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Dr Glucose, is useful for obtaining many body health data.

It is non-invasive, so you will not feel any pain during testing. You can also use it several times a day to keep track of your body health. Through countless laboratory-based experiments and combination of integrated sensor technology, Dr Glucose is also to collect various health data.

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Power Management/电量管理

- 1) Press power button to switch on Dr Glucose. 2) Dr Glucose will automatically shut down in 3 minutes without using.
- 3) Dr Glucose will automatically shut down in 10 seconds after testing.
- 4) Testing is not allowed during charging.
- 5) Press and hold power button for 3 seconds to switch off Dr Glucose.

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When Dr Glucose is switched on, LED will lit in blue.
无创测糖仪处于开机状态时, LED 会亮蓝灯

B. Ventilation area 散热区域:

The system can be cooled here.
让产品有更好的散热功效

C. Power Button 开机键:

When pressed, Dr Glucose is switched on.
启动无创测糖仪的按钮

D. Micro USB port 界面:

Dr Glucose can be charged here.
无创测糖仪的 USB 充电孔

E. Reset Button 重置键:

When pressed, Dr Glucose is reset to factory setting.
重无创测糖仪至原厂设定

F. Bluetooth LED signal 蓝芽指示灯:

When connected to a smart device, LED will lit in blue.
与产品做蓝芽链接时, LED 灯将亮蓝灯显示

- 1)Temperature range: 20 - 27 degrees
- 2)Humidity range: 20% - 70%
- 3)Atmospheric Pressure Range:700-1060hPa
- 4)Avoid direct contact with sunlight
- 5)Avoid contact with toxic or flammable gas
- 6)Working Rated Voltage: 5V
- 7)Working Rated Current: 1A

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
Please be aware of the following, or test result will be affected

- 1) Use Dr Glucose only on a flat surface.
- 2) Keep your nail up during the testing process.
- 3) Suitable finger temperature: 29 - 35 degrees
- 4) Test Interval is 10 minutes, minimum time required for heat to be fully dissipated from the sensor.
- 5) Do not apply Isolation products such as nail polish and sunblock before testing.
- 6) Allow cooling down after sports for at least 1 hour.
- 7) Allow cooling down after bathing or washing for at least 30 minutes.
- 8) Test after consuming alcohol is not allowed.
- 9) **Please put in finger into sensor, don't press down while testing otherwise it may affect test result.**
- 10) **Our APP is unable to accommodate to all existing Android smart phone models, if you encounter dropping characters or unable to input text, please set your font size to medium.**

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

1. Make sure your iOS device's Bluetooth is switched on.
2. Open up Dr Glucose app, and logged into profile, on the upper right hand corner of the app, there is an option button. Once pressed, a dropdown list menu appears.
3. Select Bluetooth.

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

Splash Page APP 开启画面

4. If the device's Bluetooth is switched on, it should be able to detect Dr Glucose device, named ITON DM: xxx. Press it to connect. Once connected, it should be highlighted in blue as shown in the picture below.
 5. You should be able to run a test once the procedure is completed.
4. 请确认蓝牙已开启, 如有开启, 用户可点选无创测糖机器代号: ITOM DM: xxx, 点选将它联结至手机,
5. 点选完毕后, ITOM DM: xxx 反蓝表示选择完毕, 用户方可开始进行测试

Every time Dr Glucose APP is opened, user will be reminded that data shown at the end of the test is for health trend monitoring and referencing only.

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
Once completed, click finish on the upper right to
complete registration.

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Login Page APP 登入画面

User will be required to login
here by choosing his profile
and password here.

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By pressing "choose user",
user can select his/her own
profile from the list.

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

By pressing the rubbish icon, user can delete unneeded profile from the APP. Password is needed in order to delete the profile.

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Insert Password 输入密码

When you have inserted incorrect password, you will not be allowed to proceed. The box will be highlighted in red.

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Once user has logged in, user enters the main menu. If the user needs to log out, he/she can do so via the upper right 'more options' button.

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Shows the height of user.(in cm)
显示用户的身高值(公分)

Testing: Enters pre-testing interface
检测: 点选进入测试前的接口

History: User can check his past test results here.
历史记录: 点选查看此用户以前的测试纪录

User Guide: User notifications can be read here.
用户指南: 点选查看用户注意事项

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Self-Calibration: User can start calibration process here.
自标定: 点选使用 30 天自我校准模式

Clear self-calibration: User can clear all self-calibration setting.
清除自标定: 点选清除所有校准的数据

Pre-testing 测试前接口

Every time the user would like to test, he/she needs to check the appropriate meal conditions and medicine (multiple choice applies), and press "begin test".

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Empty

Fasting glucose: User has not eaten anything for 8 hours or more.
清晨空腹: 用户 8 小时以上未进食

1 hour

1 hour has passed since user's previous meal. Test to see Glycemic index is high or low, varies from age and metabolism rate. The value is not real glucose. 餐后一小时: 用户离上次进食已超过一小时 测试升糖指数高或低, 其数值依据年龄及代谢速度有不同变化

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
Bluetooth connecting in progress
蓝牙设备连结中

Please insert your finger
请用户伸入手指

※Please put in finger into sensor, don't press down while testing otherwise it may affect test result.

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Environmental humidity
too low or high 环境湿
度过低或过高

Skin temperature too
low or high 体表温度
过低或过高

Data collecting error
has occurred 数据采集
异常, 错误提示


System has failed to connect to device.
搜寻设备失败

10minutes interval needed between
two tests. 测试的间隔必须至少超过
10 分钟

Sensor temp error has occurred.
测试传感器温度错误


Testing in Progress...

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After user has inserted his finger into the probe, APP will automatically start the testing process, screen will start running animation, total of six groups.

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Test Done!

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
After the test has been completed, user can remove his/her finger. The data is processing for a brief moment before showing all the health data.

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Test Results screen


At result screen, all nine health data are shown: Blood glucose level, pulse, blood flow velocity, hemoglobin, oxygen saturation, body surface temp and humidity, environmental temp and humidity. Hemoglobin is highlighted in red because it is outside normal range. Battery life is also shown.

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Via either test results screen or main screen to access the 'history' page. User may access his/her past test results to better monitor his/her health trends.

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
History Interface (Empty)

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When user has not run any test in a new day, and no data has logged yet, the bubble animation screen will appear.

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
History Interface (Graph chart)

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User can toggle two types of interface, either graph chart or data list.

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User can press on each health parameter, such as glucose, pulse, hemoglobin, speed, oxygen, to view each specific health data test results.

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Date Toggling 选择日期

On the lower right, user can switch the dates (year/month/day) to recall the test results on that specific day.

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As everyone has different body metabolism rate, there is a small part of population which the test results are largely deviated from the actual results. Their results can be self-calibrated to improve the similarity.

Step1: Choose any meal status timing (eg. fasting glucose for 8hours or more)

Step2: Measure using the user's most trusted glucose device, and insert the test result into Dr Glucose after the user has completed testing.

Step3: User must achieve 20 tests in 30 days in order to complete self-calibration process.

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Pre-Calibration Interface 校准前界面

Pre-Calibration Interface is very similar to pretesting phase, except two additions. After the fields are completed, proceed to test as per normal.

Red box A:

Test counter reminds the user to complete 20 tests in 30 days.

Red box B:

For user to insert the result obtained from user's most trusted glucose device. (in mmol)

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- 1) Before use, calibrate Dr Glucose to ensure it works normally.
- 2) Protect Dr Glucose from water. Once splattered, cut off the power supply, switch off and contact service.
- 3) Keep Dr Glucose clean. Use a soft damp cloth to clean and sterilise. DO NOT use steam to sterilise. Protect nonwaterproof finger clip probe from any liquids.
- 4) DO NOT sterilise under high temperature or disinfect with electron beam or radiation.
- 5) Ensure Dr Glucose is not disturbed by strong EMI, such as from a radio transmitter or mobile phone.
- 6) Before use, check any damage that may harm personal safety or affect Dr Glucose performance. Once per month or more frequent checks are recommended. If any damage is found, replace the damaged part immediately.
- 7) DO NOT store Dr Glucose near metal objects to avoid accidental short circuiting of the inbuilt battery.
- 8) DO NOT use Dr Glucose immediately after taking it from a cold to a warmer and more humid place.
- 9) Keep Dr Glucose from vibrating, corrosive, and flammable materials. Protect from extreme temperatures and humidity.
- 10) Dr Glucose must be maintained by authorised and qualified engineers. DO NOT disassemble.
- 11) Do not throw, hit, or scratch surface of Dr Glucose.
- 12) The power wire needs to be connected to a two-wire socket. DO NOT use a defective socket. DO NOT replace the power wire arbitrarily. Place wires, conduits, and all cables properly.
- 13) Install Dr Glucose at a fixed place to avoid injuries and damage caused by falling.
- 14) DO NOT block the Dr Glucose's ventilation to ensure effective heat dissipation.
- 15) Before testing, ensure Dr Glucose is working normally, in a suitable environment.
- 16) Check and replace frequently used accessories and cables; Dispose damaged parts as electronic waste properly.

- 17) DONOT use in a place with flammable gas, such as oxygen-mixed anesthetic.
- 18) DONOT take out the battery. DONOT throw the battery into the fire, it may explode.

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